

Blood Pressure RT (Readings and Tracking) Installation Guide


1. Click on the **Download Blood Pressure RT** button from the website.




2. Save **BloodPressureRTInstallation.zip** to your desktop.



3. Next, Right-click on BloodPressureRTInstallation.zip and select **Extract Here**.




4. Locate and open **Blood Pressure RT Installation** folder on your desktop.




5. Inside Blood Pressure RT Installation folder double-click on **setup.exe**




6. The installation wizard will be displayed. Click **Next**.




7. Read the License Agreement carefully. Select **I agree**. Click **Next** to continue.




8. Click **Next**.



9. Click **Next**.



10. Click **Close** when install is finish.



11. The installation is complete.

12. Locate on your desktop **Blood Pressure RT** icon.



13. Double-click on icon to begin.